





1.12 Data Analysis Overview

Std. Familiarization Time: 0:15:00
Approximate Video Time: 0:12:00

[Play the Video](#)

MassHunter Basics for Agilent ICP-MS
Rev.3.0.2 – February 2021
54

Batch-at-a-Glance



Batch-at-a-Glance

- [Process Batch]
- [Specify Process Order]
- [Conc] and [Count] modes
- Showing replicate data – [Detail]
- [User Columns Settings]
- [Number Configuration]
- [FullQuant Outliers]

Spectrum Pane

- ISTD Stability Graph
- QC Sample Stability Graph
- Spectrum/Chromatogram
- Adjusting the scale – right-click and drag the mouse
- Shifting the axes – left-click and drag the mouse
- Expanding a selected area – right-click and drag the mouse

Calibration Curve Pane

- Single Calibration Curve
- Twelve Calibration Curve
- [Next Element] and [Previous Element]

Method Editor

