
A video thumbnail for the '1.4 Performance Report'. It features a blue background with a white starburst logo on the left. The title '1.4 Performance Report' is in white text. Below the title, it lists 'Std. Familiarization Time: 0:05:00' and 'Approximate Video Time: 0:02:30'. A red button with white text says 'Play the Video'. At the bottom right, small text reads 'MassHunter Basics for Agilent ICP-MS Rev.3.0.2 – February 2021 16'.

## Performance Report

- Normally, the Performance Report runs at Startup
- You can monitor the status of the device status successively by executing this every time using the same conditions
- Ribbon > [Report] tab > [Performance Report]

## Performance Report



## Others

- [History View]
- [Select Data Displayed in History View]
- [Export] (Context Menu)
- [Clear Performance Report History]